



## DK50 PLUS

For one dental unit  
Pour une unité dentaire  
Para una unidad dental



Accessories/ Accessoires/ Accesorios: **BOX MD KJF FS REG AOK**

	DK50 PLUS	DK50 PLUS/M	DK50 PLUS S	DK50 PLUS S/M
*	1	1	1	1
DRYER	-	✓	-	✓
BOX	-	-	✓	✓
(5-7bar)	75 L.min <sup>-1</sup> / 5bar (50Hz) 85 L.min <sup>-1</sup> / 5bar (60Hz)	58 L.min <sup>-1</sup> / 5bar (50Hz) 68 L.min <sup>-1</sup> / 5bar (60Hz)	75 L.min <sup>-1</sup> / 5bar (50Hz) 85 L.min <sup>-1</sup> / 5bar (60Hz)	58 L.min <sup>-1</sup> / 5bar (50Hz) 68 L.min <sup>-1</sup> / 5bar (60Hz)
(6-8bar)	70 L.min <sup>-1</sup> / 6bar (50Hz) 80 L.min <sup>-1</sup> / 6bar (60Hz)	60 L.min <sup>-1</sup> / 6bar (50Hz) 70 L.min <sup>-1</sup> / 6bar (60Hz)	70 L.min <sup>-1</sup> / 6bar (50Hz) 80 L.min <sup>-1</sup> / 6bar (60Hz)	60 L.min <sup>-1</sup> / 6bar (50Hz) 70 L.min <sup>-1</sup> / 6bar (60Hz)
(8-10bar)	60 L.min <sup>-1</sup> / 8bar (50Hz) 70 L.min <sup>-1</sup> / 8bar (60Hz)	50 L.min <sup>-1</sup> / 8bar (50Hz) 60 L.min <sup>-1</sup> / 8bar (60Hz)	60 L.min <sup>-1</sup> / 8bar (50Hz) 70 L.min <sup>-1</sup> / 8bar (60Hz)	50 L.min <sup>-1</sup> / 8bar (50Hz) 60 L.min <sup>-1</sup> / 8bar (60Hz)
S1-S3	S1 - 100%	S1 - 100%	S1 - 100%	S1 - 100%
	≤ 66 dB(A)	≤ 66 dB(A)	≤ 47 dB(A)	≤ 47 dB(A)

\* Number of dental units is subject to regional requirements.

\* Le nombre d'unités dentaires dépend des exigences régionales.

\* El número de unidades dentales está sujeto a los requisitos regionales.

## DK50 2V

For two dental units  
Pour deux unités dentaires  
Para dos unidades dentales



Accessories/ Accessoires/ Accesorios: **BOX MD KJF FS REG AOK**

	DK50 2V	DK50 2V/M	DK50 2VS	DK50 2VS/M
*	2	2	2	2
DRYER	-	✓	-	✓
BOX	-	-	✓	✓
(5-7bar)	140 L.min <sup>-1</sup> / 5bar	108 L.min <sup>-1</sup> / 5bar	140 L.min <sup>-1</sup> / 5bar	108 L.min <sup>-1</sup> / 5bar
(6-8bar)	135 L.min <sup>-1</sup> / 6bar	104 L.min <sup>-1</sup> / 6bar	135 L.min <sup>-1</sup> / 6bar	104 L.min <sup>-1</sup> / 6bar
(8-10bar)	110 L.min <sup>-1</sup> / 8bar	87 L.min <sup>-1</sup> / 8bar	110 L.min <sup>-1</sup> / 8bar	87 L.min <sup>-1</sup> / 8bar
S1-S3	S1 - 100%	S1 - 100%	S1 - 100%	S1 - 100%
	≤ 67 dB(A)	≤ 69 dB(A)	≤ 51,5 dB(A)	≤ 55 dB(A)

\* Number of dental units is subject to regional requirements.

\* Le nombre d'unités dentaires dépend des exigences régionales.

\* El número de unidades dentales está sujeto a los requisitos regionales.